



### Specificaties - EL20 - Education Line pH

Korte beschrijving	Benchtop meter for pH and mV/ORP measurements
mV-bereik	-1999 to 1999
Relatieve mV-nauwkeurigheid	$\pm 1$
mV-resolutie	1
pH-bereik	0.00 to 14.00
Relatieve pH-nauwkeurigheid	$\pm 0.01$
pH-resolutie	0.01
Voeding	Ext. power supply (standard)
Temperatuurnauwkeurigheid °C	$\pm 0.5$
Temperatuurbereik °C	0.0 to 100.0
Temperatuurresolutie °C	0.1
Display	LCD
Connection to computer	Analog output
Calibration	Max. 2 points, 3 predefined buffer groups