


Wireless Datalogging for iPad® and Tablets

The NEW PASCO SPARKlink® Air wireless interface connects to desktop, laptop, and tablet computers, using either USB or Bluetooth.™




For iPad

SPARKvue® HD data acquisition and analysis software is an excellent solution for the physics lab when used with the SPARKlink Air.

PASCO sparklink® Air

- ▶ Two PASPORT™ sensor ports
- ▶ Built-in Temperature and Voltage sensors with probes
- ▶ Rechargeable lithium polymer battery
- ▶ Connect with Bluetooth or USB




For Windows® 8 tablets

Use the SPARKlink Air to connect to tablets running PASCO Capstone® software.

Shown here—Match Graph

The Motion Sensor signal reflects off the back of the tablet as the student tries to match her motion to the provided graph. The data is sent via Bluetooth from the SPARKlink Air and displayed in real time.



Order Information:

SPARKlink Air	PS-2011
PASCO Capstone	
Single User License	UI-5401
<i>For use with USB and Bluetooth on a desktop or laptop computer running Windows or Mac®, or on tablets running Windows 8</i>	
SPARKvue HD	www.pasco.com
<i>For use with Bluetooth on iPads and select Android Tablets</i>	



Store and charge up to five SPARKlink Air interfaces with a single power connection.



SPARKlink Air Charging Station PS-2577



Learn more at www.pasco.com or +1-916-786-3800